

## Datasheet

### NOL3 polyclonal antibody

**Catalog Number:** PAB13079

**Regulatory Status:** For research use only (RUO)

**Product Description:** Rabbit polyclonal antibody raised against synthetic peptide of NOL3.

**Immunogen:** A synthetic peptide corresponding to amino acids 191-208 of human NOL3.

**Host:** Rabbit

**Reactivity:** Human

**Applications:** IHC-P, WB-Ce

(See our web site product page for detailed applications information)

**Protocols:** See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Specificity:** An approximately 25 KDa band can be detected.

**Form:** Liquid

**Recommend Usage:** Western Blot (1:1000)

The optimal working dilution should be determined by the end user.

**Storage Buffer:** In PBS (0.02% sodium azide)

**Storage Instruction:** Store at 4°C for three months. For long term storage store at -20°C.

Aliquot to avoid repeated freezing and thawing.

**Entrez GeneID:** 8996

**Gene Symbol:** NOL3

**Gene Alias:** ARC, CARD2, MYC, MYP, NOP, NOP30

#### References:

1. ARC, an inhibitor of apoptosis expressed in skeletal muscle and heart that interacts selectively with caspases. Koseki T, Inohara N, Chen S, Nunez G. Proc Natl Acad Sci U S A. 1998 Apr 28;95(9):5156-60.