

Datasheet

c-Myc tag polyclonal antibody

Catalog Number: PAB13741

Regulatory Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against synthetic peptides of c-Myc tag.

Immunogen: A synthetic peptide (conjugated with KLH) corresponding to c-Myc tag.

Sequence: EQKLISEEDL

Host: Rabbit

Applications: ELISA, IF, WB

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Recommend Usage: ELISA (1:200-1:500)

Western Blot (1:2000-1:10000)

Immunofluorescence (1:250-1:1000)

The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS, pH 7.2 (0.09% sodium azide)

Storage Instruction: Store at 4°C.