

## Datasheet

### **Chlamydia trachomatis (EB) polyclonal antibody (FITC)**

**Catalog Number:** PAB14210

**Regulatory Status:** For research use only (RUO)

**Product Description:** Rabbit polyclonal antibody raised against *Chlamydia trachomatis*, Elementary Bodies.

**Immunogen:** L2 and other serovar groups.

**Host:** Rabbit

**Reactivity:** Bacteria

**Applications:** ELISA, IF

(See our web site product page for detailed applications information)

**Protocols:** See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Specificity:** Purified elementary bodies, disrupted.

Cross-reacts with *Chlamydia psittaci* and *Chlamydia pneumoniae* (TWAR). Negative vs. HEp-2 cells and egg yolk sac.

**Form:** Liquid

**Conjugation:** FITC

**Recommend Usage:** The optimal working dilution should be determined by the end user.

**Storage Buffer:** In 10 mM PBS, pH 7.2 (10 mg/mL BSA, 0.09% sodium azide)

**Storage Instruction:** Store in the dark at 4°C. For long term storage store at -20°C.

Avoid prolonged exposure to light.

Aliquot to avoid repeated freezing and thawing.