

Datasheet

Chlamydia trachomatis (EB) polyclonal antibody (FITC)

Catalog Number: PAB14211

Regulatory Status: For research use only (RUO)

Product Description: Goat polyclonal antibody raised against *Chlamydia trachomatis*, Elementary Bodies.

Immunogen: L2 and other serovar groups.

Host: Goat

Reactivity: Bacteria

Applications: ELISA, IF
(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: Purified elementary bodies, disrupted.
Cross-reacts with *Chlamydia psittacii* & *Chlamydia pneumoniae* (TWAR). Negative against HEp-2 cells and egg yolk sac. All antigens of strain L2.

Form: Liquid

Conjugation: FITC

Recommend Usage: The optimal working dilution should be determined by the end user.

Storage Buffer: In 10 mM PBS, pH 7.2 (10 mg/mL BSA, 0.09% sodium azide)

Storage Instruction: Store in the dark at 4°C. For long term storage store at -20°C.
Avoid prolonged exposure to light.
Aliquot to avoid repeated freezing and thawing.