

## Datasheet

### Diphtheria Toxin/Toxoid polyclonal antibody (FITC)

**Catalog Number:** PAB14229

**Regulatory Status:** For research use only (RUO)

**Product Description:** Goat polyclonal antibody raised against diphtheria toxoid.

**Immunogen:** Purified diphtheria toxoid.

**Host:** Goat

**Reactivity:** Bacteria

**Applications:** ELISA, GSA

(See our web site product page for detailed applications information)

**Protocols:** See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Specificity:** Demonstrated reactivity to purified diphtheria toxin/toxoid.

**Form:** Liquid

**Conjugation:** FITC

**Recommend Usage:** The optimal working dilution should be determined by the end user.

**Storage Buffer:** In 10 mM PBS, pH 7.2 (10 mg/mL BSA, 0.09% sodium azide)

**Storage Instruction:** Store in the dark at 4°C. For long term storage store at -20°C.

Avoid prolonged exposure to light.

Aliquot to avoid repeated freezing and thawing.