

Datasheet

***Escherichia coli* O and K antigens polyclonal antibody (FITC)**

Catalog Number: PAB14232

Regulatory Status: For research use only (RUO)

Product Description: Goat polyclonal antibody raised against *Escherichia coli* O and K antigens.

Immunogen: Mixture of *Escherichia coli* serotypes.

Host: Goat

Reactivity: *Escherichia coli*

Applications: ELISA, IF

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: Many "O" and "K" antigenic serotypes of *Escherichia coli*. Antiserum is not absorbed and may react with related Enterobacteriaceae.

Form: Liquid

Conjugation: FITC

Recommend Usage: The optimal working dilution should be determined by the end user.

Storage Buffer: In 10 mM PBS, pH 7.2 (10 mg/mL BSA, 0.09% sodium azide)

Storage Instruction: Store in the dark at 4°C. For long term storage store at -20°C.

Avoid prolonged exposure to light.

Aliquot to avoid repeated freezing and thawing.