

## Datasheet

### Hepatitis B virus surface antigen polyclonal antibody (FITC)

**Catalog Number:** PAB14252

**Regulatory Status:** For research use only (RUO)

**Product Description:** Goat polyclonal antibody raised against hepatitis B virus surface antigen (HBsAg) (ad/ay).

**Immunogen:** Mixture of subtypes ad & ay. hepatitis B virus surface antigen purified from human serum.

**Host:** Goat

**Reactivity:** Viruses

**Applications:** ELISA

(See our web site product page for detailed applications information)

**Protocols:** See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Specificity:** Purified surface antigen. Monospecific, reacts only with Hepatitis B virus.

**Form:** Liquid

**Conjugation:** FITC

**Recommend Usage:** The optimal working dilution should be determined by the end user.

**Storage Buffer:** In 10 mM PBS, pH 7.2 (10 mg/mL BSA, 0.09% sodium azide)

**Storage Instruction:** Store in the dark at 4°C. For long term storage store at -20°C. Avoid prolonged exposure to light. Aliquot to avoid repeated freezing and thawing.