

Datasheet

Hepatitis B virus surface antigen polyclonal antibody (FITC)

Catalog Number: PAB14254

Regulatory Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against hepatitis B virus surface antigen (HBsAg) (ad/ay).

Immunogen: Mixture of subtypes ad & ay. hepatitis B virus surface antigen purified from human serum.

Host: Rabbit

Reactivity: Viruses

Applications: IF

(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: Monospecific, reacts only with Hepatitis B surface antigen including the pre-S1 epitope.
Non-reactive with normal human serum.

Form: Liquid

Conjugation: FITC

Recommend Usage: The optimal working dilution should be determined by the end user.

Storage Buffer: In 10 mM PBS, pH 7.2 (10 mg/mL BSA, 0.09% sodium azide)

Storage Instruction: Store in the dark at 4°C. For long term storage store at -20°C.
Avoid prolonged exposure to light.
Aliquot to avoid repeated freezing and thawing.