

Datasheet

Human immunodeficiency virus type 1 p24 polyclonal antibody (FITC)

Catalog Number: PAB14271

Regulatory Status: For research use only (RUO)

Product Description: Goat polyclonal antibody raised against human immunodeficiency virus type 1 p24.

Immunogen: Native purified native p24 from strain IIIB.

Host: Goat

Reactivity: Viruses

Applications: ELISA, IF, WB

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: p24 (predominate p24 band pattern in Western blot). Antiserum does not cross-react with human T or B cells. Antiserum is not reactive with human or bovine sera.

Form: Liquid

Conjugation: FITC

Recommend Usage: The optimal working dilution should be determined by the end user.

Storage Buffer: In 10 mM PBS, pH 7.2 (10 mg/mL BSA, 0.09% sodium azide)

Storage Instruction: Store in the dark at 4°C. For long term storage store at -20°C.

Avoid prolonged exposure to light.

Aliquot to avoid repeated freezing and thawing.