

Datasheet

Influenza A (H1N1) polyclonal antibody (FITC)

Catalog Number: PAB14277

Regulatory Status: For research use only (RUO)

Product Description: Goat polyclonal antibody raised against Influenza A Strain A/USSR/90/77 (H1N1)..

Immunogen: Influenza A Strain A/USSR/90/77 (H1N1).

Host: Goat

Reactivity: Viruses

Applications: IF, IHC-P

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: Purified virions. Specific to H1N1 by IHA.

Does not react with Influenza B, RSV, Para 1-3 or Adeno. Does not react with HEp-2 cells. May react with chicken cellular proteins.

Form: Liquid

Conjugation: FITC

Recommend Usage: Recommend dilution (1:10-1:50)

The optimal working dilution should be determined by the end user.

Storage Buffer: In 10 mM PBS, pH 7.2 (10 mg/mL BSA, 0.09% sodium azide)

Storage Instruction: Store in the dark at 4°C. For long term storage store at -20°C.

Avoid prolonged exposure to light.

Aliquot to avoid repeated freezing and thawing.