

Datasheet

S-tag polyclonal antibody (FITC)

Catalog Number: PAB14866

Regulatory Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against synthetic peptide of S-tag.

Immunogen: A synthetic peptide (conjugated with FITC) corresponding to pancreatic ribonuclease A.

Sequence: KETAAAKFERQHMS

Host: Rabbit

Applications: ICC

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: This antibody recognizes protein that contain the S-tag epitope tag.

Form: Liquid

Conjugation: FITC

Recommend Usage: The optimal working dilution should be determined by the end user.

Storage Buffer: In 50 mM HEPES, 100 mM NaCl, pH 7.0 (0.2% BSA, 0.01% thimerosal)

Storage Instruction: Store at 4°C. Do not freeze.