

Datasheet

Slc7a11 polyclonal antibody (DyLight 549)

Gene Alias: 9930009M05Rik, AI451155, sut, xCT

Catalog Number: PAB14939

Regulatory Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against synthetic peptide of Slc7a11.

Immunogen: A synthetic peptide corresponding to amino acids 1-50 of mouse Slc7a11.

Host: Rabbit

Reactivity: Human

Applications: WB

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: This antibody reacts with the human protein. It does not react with murine stomach tissue in Western blot. Other species have not been tested. This is specific to XCT. By Western blot a band at ~55 KDa is seen. No other applications have been tested.

Recommended Starting Dilutions: Western Blot: 1 ug/ml (ECL) Immunoprecipitation: Not tested

Immunohisto/cytochemistry: Not tested *Optimal working dilutions should be determined by the investigator.

Form: Liquid

Conjugation: DyLight 549

Recommend Usage: The optimal working dilution should be determined by the end user.

Storage Buffer: In 50 mM sodium borate

Storage Instruction: Store at 4°C. Do not freeze.

Entrez GeneID: 26570

Gene Symbol: Slc7a11