

Datasheet

Kdr polyclonal antibody (DyLight 488)

Gene Alias: 6130401C07, Flk-1, Flk1, Krd-1, Ly73, VEGFR-2, VEGFR2

Catalog Number: PAB14948

Regulatory Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against synthetic peptide of Kdr.

Immunogen: A synthetic peptide corresponding to amino acids 1200-1300 of mouse Kdr.

Host: Rabbit

Reactivity: Human

Applications: WB

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: This antibody reacts with the human protein. Other species have not yet been tested. The immunogen used for this antibody production has 100% homology with mouse, rat and human proteins. This is specific to VEGFR-2. This antibody is useful for Western blot analysis, where a doublet is seen at ~150 KDa. Optimal working dilutions should be determined by the investigator.

Form: Liquid

Conjugation: DyLight 488

Recommend Usage: Western Blot (1:250-1:1000)

The optimal working dilution should be determined by the end user.

Storage Buffer: In 50 mM sodium borate

Storage Instruction: Store at 4°C. Do not freeze.

Entrez GeneID: 16542

Gene Symbol: Kdr