

Datasheet

SYCP1 polyclonal antibody (DyLight 488)

Gene Alias: HOM-TES-14, MGC104417, SCP1

Catalog Number: PAB14960

Regulatory Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against synthetic peptide of SYCP1.

Immunogen: A synthetic peptide corresponding to C-terminus of human SYCP1.

Host: Rabbit

Reactivity: Human, Mouse, Rat

Applications: IF

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: This SCP1 antibody recognizes human, rat and murine protein. Other species have not been tested. This is specific to human SCP1. This antibody is useful for immunofluorescence. It does not appear to work in immunohistochemistry (paraffin-embedded sections). The investigator should determine the optimal working dilution for a specific application.

Form: Liquid

Conjugation: DyLight 488

Recommend Usage: Immunofluorescence

(1:300-1:750)

The optimal working dilution should be determined by the end user.

Storage Buffer: In 50 mM sodium borate

Storage Instruction: Store at 4°C. Do not freeze.

Entrez GeneID: 6847

Gene Symbol: SYCP1