

Datasheet

Sycp1 polyclonal antibody (DyLight 549)

Catalog Number: PAB14964

Regulatory Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against synthetic peptide of Sycp1.

Immunogen: A synthetic peptide corresponding to C-terminus of mouse Sycp1.

Host: Rabbit

Reactivity: Mouse

Applications: IF, IHC-Fr, IHC-P
(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: This is specific to mouse SCP1. This antibody is useful for immunofluorescence and immunohistochemistry (paraffin and frozen sections).

Form: Liquid

Conjugation: DyLight 549

Recommend Usage: Immunofluorescence (1:100-1:750)
Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) (1:750)
Immunohistochemistry (Frozen sections) (1:750)
The optimal working dilution should be determined by the end user.

Storage Buffer: In 50 mM sodium borate

Storage Instruction: Store at 4°C. Do not freeze.

Entrez GeneID: 20957

Gene Symbol: Sycp1

Gene Alias: SCP1