

Datasheet

CYTSB polyclonal antibody (DyLight 549)

Gene Symbol: CYTSB

Gene Alias: FLJ36955, HCMOGT-1, NSP, SPECC1

Catalog Number: PAB14982

Regulatory Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against synthetic peptide of CYTSB.

Immunogen: A synthetic peptide corresponding to amino acids 200-300 of human CYTSB.

Host: Rabbit

Reactivity: Human

Applications: WB

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: This antibody reacts with human protein. There is 93% homology of the immunogen with the mouse protein and 87% homology of the immunogen with the rat NSP5a3a protein. This is specific to NSP5a3a. This antibody is useful for Western blot analysis, where a band is seen ~80 KDa. A second band may be seen ~60 KDa. The specificity of this band is not known. No other applications have been tested. Optimal working dilutions should be determined by the investigator.

Form: Liquid

Conjugation: DyLight 549

Recommend Usage: Western Blot (1:500-1:2000)

The optimal working dilution should be determined by the end user.

Storage Buffer: In 50 mM sodium borate

Storage Instruction: Store at 4°C. Do not freeze.

Entrez GeneID: 92521