

Datasheet

Lyve1 polyclonal antibody (DyLight 549)

Catalog Number: PAB14998

Regulatory Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against synthetic peptide of Lyve1.

Immunogen: A synthetic peptide corresponding to amino acids 250-318 of mouse Lyve1.

Host: Rabbit

Reactivity: Human, Mouse

Applications: Flow Cyt, IHC, WB
(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: This antibody cross-reacts with human and mouse. No other species have been tested. Rabbit Polyclonal to LYVE1 Western blot where a band at ~45 KDa is seen.

Form: Liquid

Conjugation: DyLight 549

Recommend Usage: Flow Cytometry (1:25)
The optimal working dilution should be determined by the end user.

Storage Buffer: In 50 mM sodium borate

Storage Instruction: Store at 4°C. Do not freeze.

Entrez GeneID: 114332

Gene Symbol: Lyve1

Gene Alias: 1200012G08Rik, Crsbp-1, Lyve-1, Xlkd1