

## Datasheet

### Lyve1 polyclonal antibody (DyLight 488)

**Catalog Number:** PAB14999

**Regulatory Status:** For research use only (RUO)

**Product Description:** Rabbit polyclonal antibody raised against synthetic peptide of Lyve1.

**Immunogen:** A synthetic peptide corresponding to amino acids 250-318 of mouse Lyve1.

**Host:** Rabbit

**Reactivity:** Human, Mouse

**Applications:** Flow Cyt, IHC, WB  
(See our web site product page for detailed applications information)

**Protocols:** See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Specificity:** This antibody cross-reacts with human and mouse. No other species have been tested. Rabbit Polyclonal to LYVE1 Western blot where a band at ~45 KDa is seen.

**Form:** Liquid

**Conjugation:** DyLight 488

**Recommend Usage:** Flow Cytometry (1:25)  
The optimal working dilution should be determined by the end user.

**Storage Buffer:** In 50 mM sodium borate

**Storage Instruction:** Store at 4°C. Do not freeze.

**Entrez GeneID:** 114332

**Gene Symbol:** Lyve1

**Gene Alias:** 1200012G08Rik, Crsbp-1, Lyve-1, Xlkd1