

Datasheet

CTGF polyclonal antibody (DyLight 488)

Catalog Number: PAB15039

Regulatory Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against synthetic peptide of CTGF.

Immunogen: A synthetic peptide corresponding to amino acids 299-349 of human CTGF.

Host: Rabbit

Reactivity: Mouse

Applications: WB

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: This antibody reacts with murine CTGF protein. The immunogen used for this antibody production has 100% homology with human, rat, and mouse protein. This is specific to CTGF. By Western blot a band at ~38 KDa is seen. No other applications have been tested.

Form: Liquid

Conjugation: DyLight 488

Recommend Usage: Western Blot (1:1000-1:2000)

The optimal working dilution should be determined by the end user.

Storage Buffer: In 50 mM sodium borate

Storage Instruction: Store at 4°C. Do not freeze.

Entrez GeneID: 1490

Gene Symbol: CTGF

Gene Alias: CCN2, HCS24, IGFBP8, MGC102839, NOV2