

## Datasheet

### TUBB polyclonal antibody (DyLight 649)

**Catalog Number:** PAB15098

**Regulatory Status:** For research use only (RUO)

**Product Description:** Rabbit polyclonal antibody raised against synthetic peptide of TUBB.

**Immunogen:** A synthetic peptide corresponding to amino acids 1-100 of human TUBB.

**Host:** Rabbit

**Reactivity:**

Chicken, Human, Mouse, Pig, Rat, Xenopus, Zebra fish

**Applications:** WB

(See our web site product page for detailed applications information)

**Protocols:** See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Specificity:** This is specific to beta tubulin. This antibody can be used for Western Blot where a band at ~50 KDa is seen.

**Form:** Liquid

**Conjugation:** DyLight 649

**Recommend Usage:** Western Blot (1:1000-1:10000)

The optimal working dilution should be determined by the end user.

**Storage Buffer:** In 50 mM sodium borate

**Storage Instruction:** Store at 4°C. Do not freeze.

**Entrez GeneID:** 203068

**Gene Symbol:** TUBB

**Gene Alias:** M40, MGC117247, MGC16435, OK/SW-cl.56, TUBB1, TUBB5