

Datasheet

TUBB polyclonal antibody (DyLight 488)

Catalog Number: PAB15100

Regulatory Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against synthetic peptide of TUBB.

Immunogen: A synthetic peptide corresponding to amino acids 1-100 of human TUBB.

Host: Rabbit

Reactivity:

Chicken, Human, Mouse, Pig, Rat, Xenopus, Zebra fish

Applications: WB

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: This is specific to beta tubulin. This antibody can be used for Western Blot where a band at ~50 KDa is seen.

Form: Liquid

Conjugation: DyLight 488

Recommend Usage: Western Blot (1:1000-1:10000)

The optimal working dilution should be determined by the end user.

Storage Buffer: In 50 mM sodium borate

Storage Instruction: Store at 4°C. Do not freeze.

Entrez GeneID: 203068

Gene Symbol: TUBB

Gene Alias: M40, MGC117247, MGC16435, OK/SW-cl.56, TUBB1, TUBB5