

Datasheet

Trpm2 polyclonal antibody (DyLight 649)

Catalog Number: PAB15105

Regulatory Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against synthetic peptide of Trpm2.

Immunogen: A synthetic peptide corresponding to amino acids 1200-1300 of mouse Trpm2.

Host: Rabbit

Reactivity: Mouse, Rat

Applications: WB

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: This antibody is useful in Western blot, where a band is seen ~172 KDa.

Form: Liquid

Conjugation: DyLight 649

Recommend Usage: Western Blot (2 ug/mL)

The optimal working dilution should be determined by the end user.

Storage Buffer: In 50 mM sodium borate

Storage Instruction: Store at 4°C. Do not freeze.

Entrez GeneID: 28240

Gene Symbol: Trpm2

Gene Alias: 9830168K16Rik, C79133, LTRPC2, TRPC7, Trp7, Trpp7