

Datasheet

Trpm2 polyclonal antibody (DyLight 488)

Catalog Number: PAB15108

Regulatory Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against synthetic peptide of Trpm2.

Immunogen: A synthetic peptide corresponding to amino acids 1430-1508 of rat Trpm2.

Host: Rabbit

Reactivity: Mouse, Rat

Applications: IHC, WB

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: This is specific to isoform 1, or the long isoform (TRPM2-L), of TRPM2. This antibody is useful for Western blot, where a band is seen ~181 KDa.

Form: Liquid

Conjugation: DyLight 488

Recommend Usage: Immunohistochemistry (1:200)

Western Blot (1:500-1:1000)

The optimal working dilution should be determined by the end user.

Storage Buffer: In 50 mM sodium borate

Storage Instruction: Store at 4°C. Do not freeze.

Entrez GeneID: 294329

Gene Symbol: Trpm2

Gene Alias: Trpm2-predicted