

## Datasheet

### TERF2 polyclonal antibody (DyLight 649)

**Catalog Number:** PAB15118

**Regulatory Status:** For research use only (RUO)

**Product Description:** Rabbit polyclonal antibody raised against recombinant TERF2.

**Immunogen:** Recombinant protein corresponding to TERF2.

**Host:** Rabbit

**Reactivity:** Human

**Applications:** ChIP, IF, IP, WB

(See our web site product page for detailed applications information)

**Protocols:** See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Specificity:** TRF2 This antibody is useful for Western blot, where a band at ~56 KDa is seen.

**Form:** Liquid

**Conjugation:** DyLight 649

**Recommend Usage:** ChIP (1:2000-1:5000)

Immunofluorescence (1:2000-1:5000)

Western Blot (1:2000-1:5000)

The optimal working dilution should be determined by the end user.

**Storage Buffer:** In 50 mM sodium borate

**Storage Instruction:** Store at 4°C. Do not freeze.