

Datasheet

TERF2 polyclonal antibody (DyLight 488)

Catalog Number: PAB15120

Regulatory Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against recombinant TERF2.

Immunogen: Recombinant protein corresponding to TERF2.

Host: Rabbit

Reactivity: Human

Applications: ChIP, IF, IP, WB

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: TRF2 This antibody is useful for Western blot, where a band at ~56 KDa is seen.

Form: Liquid

Conjugation: DyLight 488

Recommend Usage: ChIP (1:2000-1:5000)

Immunofluorescence (1:2000-1:5000)

Western Blot (1:2000-1:5000)

The optimal working dilution should be determined by the end user.

Storage Buffer: In 50 mM sodium borate

Storage Instruction: Store at 4°C. Do not freeze.