

Datasheet

Plin3 polyclonal antibody (DyLight 649)

Catalog Number: PAB15145

Regulatory Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against synthetic peptide of Plin3.

Immunogen: A synthetic peptide corresponding to amino acids 1-100 of mouse Plin3.

Host: Rabbit

Reactivity: Human, Mouse

Applications: WB

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: TIP47 This antibody is useful for Western blot analysis where a band is seen ~47 KDa (isoform B, containing the first 184 residues that isoform A is lacking).

Form: Liquid

Conjugation: DyLight 649

Recommend Usage: Western Blot (1:1000-1:5000)

The optimal working dilution should be determined by the end user.

Storage Buffer: In 50 mM sodium borate

Storage Instruction: Store at 4°C. Do not freeze.

Entrez GeneID: 66905

Gene Symbol: M6prbp1

Gene Alias: 1300012C15Rik, Tip47