

Datasheet

SPDYA polyclonal antibody (DyLight 649)

SPY1

Catalog Number: PAB15206**Regulatory Status:** For research use only (RUO)**Product Description:** Rabbit polyclonal antibody raised against synthetic peptide of SPDYA.**Immunogen:** A synthetic peptide corresponding to amino acids 200-286 of human SPDYA.**Host:** Rabbit**Reactivity:** Mouse**Applications:** WB

(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols**Specificity:** This is specific to Speedy/Ringo protein. This antibody is useful for Western blot analysis, where a band is seen at ~36 KDa (isoform 2). After a longer exposure, a 31 KDa band can be detected (isoform 1). There is also a non-specific band at ~58 KDa. This antibody was tested against human and mouse brain samples, but failed to detect proteins of the predicted size in either lysate.**Form:** Liquid**Conjugation:** DyLight 649**Recommend Usage:** The optimal working dilution should be determined by the end user.**Storage Buffer:** In 50 mM sodium borate**Storage Instruction:** Store at 4°C. Do not freeze.**Entrez GeneID:** 245711**Gene Symbol:** SPDYA**Gene Alias:** MGC110856, MGC57218, Ringo3, SPDY1,