

Datasheet

SLIRP polyclonal antibody (DyLight 549)

Catalog Number: PAB15212

Regulatory Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against synthetic peptide of SLIRP.

Immunogen: A synthetic peptide corresponding to C-terminus of human SLIRP.

Host: Rabbit

Reactivity: Human, Primates

Applications: WB

(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: Does not react with mouse. This is specific to SLIRP protein. This antibody may be used for Western blot analysis where a band at ~12 KDa is seen, representing the SLIRP protein. May see additional bands above 45 KDa.

Form: Liquid

Conjugation: DyLight 549

Recommend Usage: Western Blot (0.5-2 ug/mL)
The optimal working dilution should be determined by the end user.

Storage Buffer: In 50 mM sodium borate

Storage Instruction: Store at 4°C. Do not freeze.

Entrez GeneID: 81892

Gene Symbol: SLIRP

Gene Alias: DC23, C14orf156, DC50, PD04872

Gene Summary: Steroid receptor RNA activator (SRA,

or SRA1; MIM 603819) is a complex RNA molecule containing multiple stable stem-loop structures that functions in coactivation of nuclear receptors. SLIRP interacts with stem-loop structure-7 of SRA (STR7) and modulates nuclear receptor transactivation (Hatchell et al., 2006 [PubMed 16762838]).[supplied by OMIM]