

Datasheet

SYCP3 polyclonal antibody (DyLight 649)

Catalog Number: PAB15238

Regulatory Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against synthetic peptide of SYCP3.

Immunogen: A synthetic peptide corresponding to C-terminus of human SYCP3.

Host: Rabbit

Reactivity: Human, Mouse

Applications: IF, IHC-P, WB

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: This SCP3 antibody recognizes human and murine protein. Other species have not been tested.

Form: Liquid

Conjugation: DyLight 649

Recommend Usage: Western Blot (0.5 ug/mL)

The optimal working dilution should be determined by the end user.

Storage Buffer: In 50 mM sodium borate

Storage Instruction: Store at 4°C. Do not freeze.

Entrez GeneID: 50511

Gene Symbol: SYCP3

Gene Alias: COR1, MGC71888, SCP3