

Datasheet

Sycp3 polyclonal antibody (DyLight 649)

Catalog Number: PAB15245

Regulatory Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against synthetic peptide of Sycp3.

Immunogen: A synthetic peptide corresponding to C-terminus of mouse Sycp3.

Host: Rabbit

Reactivity: Mouse

Applications: IF, IHC-P

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: This SCP3 antibody recognizes mouse protein. Other species have not been tested. This is specific to mouse SCP3.

Form: Liquid

Conjugation: DyLight 649

Recommend Usage: Immunofluorescence

(1:100-1:750)

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) (1:750)

The optimal working dilution should be determined by the end user.

Storage Buffer: In 50 mM sodium borate

Storage Instruction: Store at 4°C. Do not freeze.

Entrez GeneID: 20962

Gene Symbol: Sycp3

Gene Alias: Cor1, Scp3