

## Datasheet

### Sycp3 polyclonal antibody (DyLight 549)

**Catalog Number:** PAB15246

**Regulatory Status:** For research use only (RUO)

**Product Description:** Rabbit polyclonal antibody raised against synthetic peptide of Sycp3.

**Immunogen:** A synthetic peptide corresponding to C-terminus of mouse Sycp3.

**Host:** Rabbit

**Reactivity:** Mouse

**Applications:** IF, IHC-P

(See our web site product page for detailed applications information)

**Protocols:** See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Specificity:** This SCP3 antibody recognizes mouse protein. Other species have not been tested. This is specific to mouse SCP3.

**Form:** Liquid

**Conjugation:** DyLight 549

**Recommend Usage:** Immunofluorescence (1:100-1:750)

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) (1:750)

The optimal working dilution should be determined by the end user.

**Storage Buffer:** In 50 mM sodium borate

**Storage Instruction:** Store at 4°C. Do not freeze.

**Entrez GeneID:** 20962

**Gene Symbol:** Sycp3

**Gene Alias:** Cor1, Scp3