

Datasheet

Trpv1 polyclonal antibody

Catalog Number: PAB16889

Regulatory Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against synthetic peptide of Trpv1.

Immunogen: A synthetic peptide (conjugated with KLH) corresponding to Trpv1.

Sequence: GEKPPRLYDRR

Host: Rabbit

Reactivity: Mouse

Applications: EIA, WB

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Lyophilized

Recommend Usage: EIA (1:100000)

Western Blot (1:5000)

The optimal working dilution should be determined by the end user.

Storage Buffer: Lyophilized from 10 mM PBS

Storage Instruction: Store at -20°C on dry atmosphere.

After reconstitution with deionized water, store at -20°C or lower.

Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 193034

Gene Symbol: Trpv1

Gene Alias: OTRPC1, TRPV1alpha, TRPV1beta, VR-1, Vr1