

Datasheet

c-Myc tag polyclonal antibody

Catalog Number: PAB17055

Regulatory Status: For research use only (RUO)

Product Description: Goat polyclonal antibody raised against synthetic peptides of c-Myc tag.

Immunogen: A synthetic peptide (conjugated with KLH) corresponding to c-Myc tag.

Sequence: EQKLISEEDL

Host: Goat

Applications: ELISA, ICC, IHC, WB
(See our web site product page for detailed applications information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Recommend Usage: ELISA (2-10 ug/mL)
Immunohistochemistry (0.5-5 ug/mL)
Immunocytochemistry (0.5-5 ug/mL)
The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS

Storage Instruction: Store at 4°C.