

Datasheet

TRPM8 polyclonal antibody

Catalog Number: PAB18678

Regulatory Status: For research use only (RUO)

Product Description: Goat polyclonal antibody raised against synthetic peptide of TRPM8.

Immunogen: A synthetic peptide corresponding to amino acids at internal region of human TRPM8.

Sequence: C-TSEQDKDNWNGQLK

Host: Goat

Applications: ELISA

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Purification: Antigen affinity purification

Concentration: 0.5 mg/mL

Recommend Usage: ELISA (1:2000)

The optimal working dilution should be determined by the end user.

Storage Buffer: In Tris saline, pH 7.3 (0.02% sodium azide, 0.5% BSA)

Storage Instruction: Store at -20°C.

Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 79054

Gene Symbol: TRPM8

Gene Alias: LTRPC6, MGC2849, TRPP8