

## Datasheet

### THY1 polyclonal antibody

**Catalog Number:** PAB18705

**Regulatory Status:** For research use only (RUO)

**Product Description:** Goat polyclonal antibody raised against synthetic peptide of THY1.

**Immunogen:** A synthetic peptide corresponding to amino acids at internal region of human THY1.

**Sequence:** C-PEHTYRSRTNFTSKY

**Host:** Goat

**Theoretical MW (kDa):** 26

**Reactivity:** Human

**Applications:** ELISA, WB-Ti  
(See our web site product page for detailed applications information)

**Protocols:** See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Form:** Liquid

**Purification:** Antigen affinity purification

**Concentration:** 0.5 mg/mL

**Recommend Usage:** ELISA (1:64000)  
Western Blot (0.1-0.3 ug/mL)  
The optimal working dilution should be determined by the end user.

**Storage Buffer:** In Tris saline, pH 7.3 (0.02% sodium azide, 0.5% BSA)

**Storage Instruction:** Store at -20°C.  
Aliquot to avoid repeated freezing and thawing.

**Entrez GeneID:** 7070

**Gene Symbol:** THY1

**Gene Alias:** CD90, FLJ33325