

## Datasheet

### POU5F1 polyclonal antibody

**Catalog Number:** PAB19184

**Regulatory Status:** For research use only (RUO)

**Product Description:** Rabbit polyclonal antibody raised against synthetic peptide of POU5F1.

**Immunogen:** A synthetic peptide corresponding to 14 amino acids near C-terminus of human POU5F1.

**Host:** Rabbit

**Reactivity:** Human, Mouse, Rat

**Applications:** ELISA, IHC, WB-Ti  
(See our web site product page for detailed applications information)

**Protocols:** See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Form:** Liquid

**Purification:** Affinity purification

**Concentration:** 1 mg/mL

**Recommend Usage:** Western blot (1-2 ug/mL)  
The optimal working dilution should be determined by the end user.

**Storage Buffer:** In PBS (0.02% sodium azide)

**Storage Instruction:** Store at 4°C for three months. For long term storage store at -20°C.  
Aliquot to avoid repeated freezing and thawing.

**Entrez GeneID:** 5460

**Gene Symbol:** POU5F1

**Gene Alias:** MGC22487, OCT3, OCT4, OTF3, OTF4