

Datasheet

SOX2 polyclonal antibody

Catalog Number: PAB19185

Regulatory Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against synthetic peptide of SOX2.

Immunogen: A synthetic peptide corresponding to 15 amino acids near N-terminus of human SOX2.

Host: Rabbit

Reactivity: Human, Mouse, Rat

Applications: ELISA, ICC, WB-Ce

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Purification: Peptide affinity purification

Concentration: 1 mg/mL

Recommend Usage: Western Blot (1-2 ug/mL)

Immunocytochemistry (5 ug/mL)

The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS (0.02% sodium azide)

Storage Instruction: Store at 4°C for three months. For long term storage store at -20°C.

Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 6657

Gene Symbol: SOX2

Gene Alias: ANOP3, MCOPS3, MGC2413

Gene Summary: This intronless gene encodes a member of the SRY-related HMG-box (SOX) family of transcription factors involved in the regulation of

embryonic development and in the determination of cell fate. The product of this gene is required for stem-cell maintenance in the central nervous system, and also regulates gene expression in the stomach. Mutations in this gene have been associated with optic nerve hypoplasia and with syndromic microphthalmia, a severe form of structural eye malformation. This gene lies within an intron of another gene called SOX2 overlapping transcript (SOX2OT). [provided by RefSeq]