

Datasheet

NOL3 polyclonal antibody

Catalog Number: PAB19787

Regulatory Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against synthetic peptide of NOL3.

Immunogen: A synthetic peptide corresponding to 14 amino acids at N-terminus of human NOL3.

Host: Rabbit

Reactivity: Human

Applications: ELISA, WB-Ce

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Recommend Usage: ELISA (1:10000-1:40000)

Western Blot (1:200-500)

The optimal working dilution should be determined by the end user.

Storage Buffer: In serum (0.02% sodium azide)

Storage Instruction: Store at 4°C for three months. For long term storage store at -20°C.

Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 8996

Gene Symbol: NOL3

Gene Alias: ARC, CARD2, MYC, MYP, NOP, NOP30