

Datasheet

BrdU polyclonal antibody

Catalog Number: PAB19944

Regulatory Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against recombinant BrdU.

Immunogen: BromodeoxyUridine conjugated with KLH.

Host: Rabbit

Applications: IF

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Purification: Immunoaffinity purification

Recommend Usage: Immunofluorescence (1:500)

The optimal working dilution should be determined by the end user.

Storage Buffer: In 20 mM KH_2PO_4 , 150 mM NaCl, pH 7.2 (0.01% sodium azide)

Storage Instruction: Store at 4°C. For long term storage store at -20°C.

Aliquot to avoid repeated freezing and thawing.