

Datasheet

H3K27me3 polyclonal antibody

Catalog Number: PAB24549

Regulatory Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against H3K27me3.

Immunogen: A synthetic peptide (conjugated with KLH) corresponding to histone H3 containing the trimethylated lysine 27 (H3K27me3).

Host: Rabbit

Reactivity: Human

Applications: ChIP, Dot, ELISA, WB
(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Purification: Affinity purification

Recommend Usage: ELISA (1:1000)
Dot Blot (1:10000)
Western Blot (1:100)
Immunofluorescence (1:200)
ChIP (1 ug/ChIP)
The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS (0.05% sodium azide, 0.05% ProClin 300)

Storage Instruction: Store at -20°C. For long term storage store at -80°C.
Aliquot to avoid repeated freezing and thawing.