

Datasheet

ESRRB polyclonal antibody

Catalog Number: PAB25698

Regulatory Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against synthetic peptide of ESRRB.

Immunogen: A synthetic peptide corresponding to 19 amino acids at C-terminus of human ESRRB.

Host: Rabbit

Reactivity: Human

Applications: ELISA, IHC, WB-Ti
(See our web site product page for detailed applications information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Purification: Peptide affinity purification

Concentration: 1 mg/mL

Isotype: IgG

Recommend Usage: Western Blot (1 ug/mL)
Immunohistochemistry (5 ug/mL)
The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS (0.05% sodium azide)

Storage Instruction: Store at 4°C for three months. For long term storage store at -20°C.
Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 2103

Gene Symbol: ESRRB

Gene Alias: DFNB35, ERR2, ERRb, ERRbeta, ERRbeta-2, ESRL2, NR3B2

Gene Summary: This gene encodes a protein with

similarity to the estrogen receptor. Its function is unknown; however, a similar protein in mouse plays an essential role in placental development. [provided by RefSeq]