

Datasheet

NOL3 polyclonal antibody

Gene Alias: ARC, CARD2, MYC, MYP, NOP, NOP30

Catalog Number: PAB27021

Regulatory Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against synthetic peptide of NOL3.

Immunogen: A synthetic peptide corresponding to human NOL3.

Host: Rabbit

Theoretical MW (kDa): 23

Reactivity: Human, Mouse

Applications: ELISA, IHC-P, WB-Ce
(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: NOL3 polyclonal antibody detects endogenous levels of NOL3 protein.

Form: Liquid

Purification: Antigen affinity purification

Concentration: 1 mg/mL

Recommend Usage: Western Blot (1:500-1:1000)
Immunohistochemistry (1:50-1:200)
Immunofluorescence (1:50-1:200)
The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS, pH 7.2 (0.05% sodium azide)

Storage Instruction: Store at 4°C. For long term storage store at -20°C.
Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 8996

Gene Symbol: NOL3