

Datasheet

ESRRB polyclonal antibody

ERRbeta-2, ESRL2, NR3B2

Catalog Number: PAB27787**Regulatory Status:** For research use only (RUO)**Product Description:** Rabbit polyclonal antibody raised against synthetic peptide of ESRRB.**Immunogen:** A synthetic peptide corresponding to 15 amino acid at ligand-binding domain of human ESRRB.**Host:** Rabbit**Reactivity:**

Bovine,Dog,Hamster,Horse,Human,Monkey,Pig,Rabbit

Applications: IHC-P

(See our web site product page for detailed applications information)

Protocols: See our web site at<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols**Specificity:** BLAST analysis of the peptide immunogen showed no homology with other human proteins, except ESRRG (67%).**Form:** Liquid**Purification:** Immunoaffinity chromatography**Recommend Usage:** Immunohistochemistry

(Formalin/PFA-fixed paraffin-embedded sections) (3 ug/mL)

The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS (0.09% sodium azide)**Storage Instruction:** Store at 4°C. For long term storage store at -80°C.

Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 2103**Gene Symbol:** ESRRB**Gene Alias:** DFNB35, ERR2, ERRb, ERRbeta,**Gene Summary:** This gene encodes a protein with similarity to the estrogen receptor. Its function is unknown; however, a similar protein in mouse plays an essential role in placental development. [provided by RefSeq]