

## Datasheet

### TRPM8 polyclonal antibody

**Catalog Number:** PAB27834

**Regulatory Status:** For research use only (RUO)

**Product Description:** Rabbit polyclonal antibody raised against synthetic peptide of TRPM8.

**Immunogen:** A synthetic peptide corresponding to 15 amino acid at internal region of human TRPM8.

**Host:** Rabbit

**Reactivity:** Chicken, Dog, Guinea pig, Hamster, Horse, Human, Mouse, Rabbit, Rat

**Applications:** IHC-P  
(See our web site product page for detailed applications information)

**Protocols:** See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Specificity:** BLAST analysis of the peptide immunogen showed no homology with other human proteins, except TRPM2 (60%).

**Form:** Liquid

**Purification:** Immunoaffinity chromatography

**Recommend Usage:** Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) (2.4 ug/mL)  
The optimal working dilution should be determined by the end user.

**Storage Buffer:** In PBS (0.09% sodium azide)

**Storage Instruction:** Store at 4°C. For long term storage store at -80°C.  
Aliquot to avoid repeated freezing and thawing.

**Entrez GeneID:** 79054

**Gene Symbol:** TRPM8

**Gene Alias:** LTRPC6, MGC2849, TRPP8