

Datasheet

SLC12A3 polyclonal antibody (FITC)

Gene Alias: FLJ96318, NCCT, TSC

Catalog Number: PAB28760

Regulatory Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against recombinant Human, Rat SLC12A3.

Immunogen: Recombinant protein corresponding to SLC12A3 N-terminus (amino acids 74-95).

Host: Rabbit

Reactivity: Human, Mouse, Rat

Applications: IF, WB

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: 1 ug/mL of PAB28760 was sufficient for detection of SLC12A3, ~160kDa, in 10ug of rat kidney tissue lysate by colorimetric immunoblot analysis using Goat anti-rabbit IgG:HRP as the secondary antibody.

Form: Liquid

Conjugation: FITC

Purification: Affinity purification

Recommend Usage: Immunofluorescence

Western Blot (1:1000)

The optimal working dilution should be determined by the end user.

Storage Buffer: PBS, 50% glycerol and 0.09% sodium azide

Storage Instruction: Store at 4°C.

Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 6559|54300

Gene Symbol: SLC12A3