

Datasheet

Scnn1a polyclonal antibody (PerCP)

Gene Alias: -

Catalog Number: PAB28781

Regulatory Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against recombinant Rat Scnn1a.

Immunogen: Recombinant protein corresponding to Scnn1a N-terminus (amino acids 46-68).

Host: Rabbit

Reactivity: Mouse,Rat

Applications: IF, IHC, WB
(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: 1 ug/mL of PAB28781 was sufficient for detection of Scnn1a, ~85kDa, in 35ug of rat kidney tissue lysate by colorimetric immunoblot analysis using Goat anti-rabbit IgG:HRP as the secondary antibody.

Form: Liquid

Conjugation: PerCP

Purification: Affinity purification

Recommend Usage: Immunofluorescence (1:400)

Immunohistochemistry

Western Blot (1:1000)

The optimal working dilution should be determined by the end user.

Storage Buffer: PBS, 50% glycerol and 0.09% sodium azide

Storage Instruction: Store at 4°C.

Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 25122

Gene Symbol: Scnn1a