

Datasheet

Scnn1b polyclonal antibody (PerCP)

Gene Symbol: Scnn1b

Catalog Number: PAB28799

Gene Alias: RNENACB

Regulatory Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against recombinant Rat Scnn1b.

Immunogen: Recombinant protein corresponding to Scnn1b C-terminus (amino acids 617-638).

Host: Rabbit

Reactivity: Human, Mouse, Rat

Applications: IF, IHC, IP, WB
(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: 1 ug/mL of PAB28799 was sufficient for detection of Scnn1b, ~87kDa, in 20ug of rat kidney tissue lysate by colorimetric immunoblot analysis using Goat anti-rabbit IgG:HRP as the secondary antibody.

Form: Liquid

Conjugation: PerCP

Purification: Affinity purification

Recommend Usage: Immunofluorescence

Immunohistochemistry

Immunoprecipitation

Western Blot (1:1000)

The optimal working dilution should be determined by the end user.

Storage Buffer: PBS, 50% glycerol and 0.09% sodium azide

Storage Instruction: Store at 4°C.
Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 24767