

## Datasheet

### Scnn1g polyclonal antibody (PerCP)

Gene Alias: ENaC

**Catalog Number:** PAB28817

**Regulatory Status:** For research use only (RUO)

**Product Description:** Rabbit polyclonal antibody raised against recombinant Rat Scnn1g.

**Immunogen:** Recombinant protein corresponding to Scnn1g C-terminus (amino acids 629-650).

**Host:** Rabbit

**Reactivity:** Mouse, Rat

**Applications:** IF, IHC, WB

(See our web site product page for detailed applications information)

**Protocols:** See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Specificity:** 1 ug/mL of PAB28817 was sufficient for detection of Scnn1g, ~83kDa, in 20ug of rat kidney tissue lysate by colorimetric immunoblot analysis using Goat anti-rabbit IgG:HRP as the secondary antibody.

**Form:** Liquid

**Conjugation:** PerCP

**Purification:** Affinity purification

**Recommend Usage:** Immunofluorescence

Immunohistochemistry

Western Blot (1:1000)

The optimal working dilution should be determined by the end user.

**Storage Buffer:** PBS, 50% glycerol and 0.09% sodium azide

**Storage Instruction:** Store at 4°C.

Aliquot to avoid repeated freezing and thawing.

**Entrez GeneID:** 24768

**Gene Symbol:** Scnn1g