

## Datasheet

### Slc14a2 polyclonal antibody (PerCP)

Gene Alias: Slc14a2T, UrT1-C, UrT1-D, testSymbol

**Catalog Number:** PAB28835

**Regulatory Status:** For research use only (RUO)

**Product Description:** Rabbit polyclonal antibody raised against recombinant Rat Slc14a2.

**Immunogen:** Recombinant protein corresponding to Slc14a2 C-terminus (amino acids 911-929).

**Host:** Rabbit

**Reactivity:** Mouse, Rat

**Applications:** IF, WB

(See our web site product page for detailed applications information)

**Protocols:** See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Specificity:** 1 ug/mL of PAB28835 was sufficient for detection of Slc14a2, ~97 and 127kDa, in 20ug of rat kidney tissue lysate by colorimetric immunoblot analysis using Goat anti-rabbit IgG:HRP as the secondary antibody.

**Form:** Liquid

**Conjugation:** PerCP

**Purification:** Affinity purification

**Recommend Usage:** Immunofluorescence

Western Blot (1:1000)

The optimal working dilution should be determined by the end user.

**Storage Buffer:** PBS, 50% glycerol and 0.09% sodium azide

**Storage Instruction:** Store at 4°C.

Aliquot to avoid repeated freezing and thawing.

**Entrez GeneID:** 54302

**Gene Symbol:** Slc14a2