

## Datasheet

### Hsp90b1 polyclonal antibody (PerCP)

**Gene Alias:** ERp99, GRP94, TA-3, Targ2, Tra-1, Tra1,  
gp96

**Catalog Number:** PAB28961

**Regulatory Status:** For research use only (RUO)

**Product Description:** Rabbit polyclonal antibody raised against synthetic peptide of Hsp90b1.

**Immunogen:** A synthetic peptide corresponding to amino acids at C-terminus of mouse Hsp90b1.

**Host:** Rabbit

**Reactivity:** Bovine, Human, Mouse, Rat

**Applications:** ICC, IHC, IP, WB

(See our web site product page for detailed applications information)

**Protocols:** See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Specificity:** Detects ~94 kDa.

**Form:** Liquid

**Conjugation:** PerCP

**Purification:** Peptide Affinity Purified

**Recommend Usage:** Immunocytochemistry

Immunohistochemistry

Immunoprecipitation (1:80)

Western blot (1:1000)

The optimal working dilution should be determined by the end user.

**Storage Buffer:** In PBS, pH 7.4 (50% glycerol, 0.09% sodium azide)

**Storage Instruction:** Store at 4°C.

Aliquot to avoid repeated freezing and thawing.

**Entrez GeneID:** 22027

**Gene Symbol:** Hsp90b1